True-Fit DENTAL IOS

Remake Scan Quick Reference

REMAKE CASE WORKFLOW · STEP 3

PROBLEM	SCAN(S) TO TAKE
Contacts	Scan the adjacent teeth with the prep. **No cord needed.** Send only the scan needed.
Occlusion high or light	New **bite scan**, full and complete.
Margin or prep fit open or short margin	**Pack cord** to fully expose the margin. Complete **new scan** if at all practical.*
Fracture	**Pack cord** to fully expose the margin. Complete **new scan** if at all practical.*
Implant rotation, contact shift, or occlusal discrepancy	New **scan body scan**: complete quadrant upper, lower, and bite. Consider a new scan body.
Contour or length	Scan with the **crown in place**, showing the adjacent and contralateral teeth.

*A die scan with cord packed is risky but can work if necessary.

True-Fit Dental IOS | Trusted lab partner: Digital Dental Leaders
Eric True | **951-430-9888** | truefitlab@gmail.com | **truefitdentalios.com**

True-Fit DENTAL IOS

Before You Rescan

REMAKE CASE WORKFLOW · STEPS 1 & 4

Every remake: notes for the lab with the specific reason and your recommended solution, plus pictures of the problem.

PROBLEM	PICTURE TO TAKE
Open or short margin	Crown in place
Open contact	Ideally with a measurement tool visible
Shade	Crown in the mouth with the shade tab visible
Occlusion	Showing the high or light occlusion
Fracture	Showing the fracture

Send the case with:

- Improvements made
- Notes in the Rx
- Pictures of the problem
- Scan(s) from the front of this card

True-Fit Dental IOS | Trusted lab partner: Digital Dental Leaders
Eric True | **951-430-9888** | truefitlab@gmail.com | **truefitdentalios.com**