True-Fit DENTAL IOS

Remake Scan Sheet

REMAKE CASE WORKFLOW · QUICK REFERENCE

> Every remake: notes for the lab with the specific reason and your recommended solution, plus pictures of the problem.

1 PICTURES TO TAKE FIRST

PROBLEM	PICTURE TO TAKE
Open or short margin	Crown in place
Open contact	Ideally with a measurement tool visible
Shade	Crown in the mouth, shade tab visible
Occlusion	Showing the high or light spot
Fracture	Showing the fracture

2 SCAN(S) TO TAKE

PROBLEM	SCAN(S) TO TAKE
Contacts	Adjacent teeth with the prep. **No cord needed.** Send only the scan needed.
Occlusion high or light	New **bite scan**, full and complete.
Margin or prep fit open or short margin	**Pack cord** to fully expose the margin; complete **new scan** if at all practical.*
Fracture	**Pack cord** to fully expose the margin; complete **new scan** if at all practical.*
Implant rotation, contact shift or occlusal discrepancy	New **scan body scan**: complete quadrant upper, lower and bite. Consider a new scan body.
Contour or length	Scan with the **crown in place**, showing adjacent and contralateral teeth.

*A die scan with cord packed is risky but can work if necessary.

3 SEND THE CASE WITH

- Improvements made
- Notes in the Rx
- Pictures of the problem
- The scan(s) from step 2

True-Fit Dental IOS | Trusted lab partner: Digital Dental Leaders
Eric True | **951-430-9888** | truefitlab@gmail.com | **truefitdentalios.com**

True-Fit DENTAL IOS

Remake Scan Sheet

REMAKE CASE WORKFLOW · QUICK REFERENCE

> Every remake: notes for the lab with the specific reason and your recommended solution, plus pictures of the problem.

1 PICTURES TO TAKE FIRST

PROBLEM	PICTURE TO TAKE
Open or short margin	Crown in place
Open contact	Ideally with a measurement tool visible
Shade	Crown in the mouth, shade tab visible
Occlusion	Showing the high or light spot
Fracture	Showing the fracture

2 SCAN(S) TO TAKE

PROBLEM	SCAN(S) TO TAKE
Contacts	Adjacent teeth with the prep. **No cord needed.** Send only the scan needed.
Occlusion high or light	New **bite scan**, full and complete.
Margin or prep fit open or short margin	**Pack cord** to fully expose the margin; complete **new scan** if at all practical.*
Fracture	**Pack cord** to fully expose the margin; complete **new scan** if at all practical.*
Implant rotation, contact shift or occlusal discrepancy	New **scan body scan**: complete quadrant upper, lower and bite. Consider a new scan body.
Contour or length	Scan with the **crown in place**, showing adjacent and contralateral teeth.

*A die scan with cord packed is risky but can work if necessary.

3 SEND THE CASE WITH

- Improvements made
- Notes in the Rx
- Pictures of the problem
- The scan(s) from step 2

True-Fit Dental IOS | Trusted lab partner: Digital Dental Leaders
Eric True | **951-430-9888** | truefitlab@gmail.com | **truefitdentalios.com**