True-Fit DENTAL IOS

Crown & Bridge Case Checklist

CHAIRSIDE · Prepared by Eric True

Patient _____ **Chart #** _____

Date _____ **Tooth #(s)** _____ **Shade** _____

Doctor _____ **Assistant initials** _____

Check 1: Preprep

- [ ] Did we take the shade following our workflow?
- [ ] Did we take the prescan to help with the occlusal check and for future reference?

Check 2: Prescan

- [ ] Did we check the prep for adequate reduction, clear margins, and undercuts?
- [ ] Did we adjust the contacts to make them parallel, broad, and smooth?
- [ ] Did we measure the occlusal clearance using the prescan + initial prep scan?
- [ ] Did we double pack cord, follow CDC, and pull the top cord?

Check 3: Postscan

- [ ] Did we check the scan fully for clear margins, complete contacts, and correct occlusion?
- [ ] Did we add to, adjust, and polish the temporary for an ideal fit?
- [ ] Did we include all the info needed in the Rx?

Notes Dry-erase marker · wipe clean after each patient

Scan to visit
True-Fit Dental IOS

True-Fit Dental IOS | Trusted lab partner: Digital Dental Leaders
Eric True | **951-430-9888** | truefitlab@gmail.com | **truefitdentalios.com**

True-Fit DENTAL IOS

Crown & Bridge Case Checklist

CHAIRSIDE · Prepared by Eric True

Patient _____ Chart # _____

Date _____ Tooth #(s) _____ Shade _____

Doctor _____ Assistant initials _____

Check 1: Preprep

- [ ] Did we take the shade following our workflow?
- [ ] Did we take the prescan to help with the occlusal check and for future reference?

Check 2: Prescan

- [ ] Did we check the prep for adequate reduction, clear margins, and undercuts?
- [ ] Did we adjust the contacts to make them parallel, broad, and smooth?
- [ ] Did we measure the occlusal clearance using the prescan + initial prep scan?
- [ ] Did we double pack cord, follow CDC, and pull the top cord?

Check 3: Postscan

- [ ] Did we check the scan fully for clear margins, complete contacts, and correct occlusion?
- [ ] Did we add to, adjust, and polish the temporary for an ideal fit?
- [ ] Did we include all the info needed in the Rx?

Notes Dry-erase marker · wipe clean after each patient

Scan to visit
True-Fit Dental IOS

True-Fit Dental IOS | Trusted lab partner: Digital Dental Leaders
Eric True | 951-430-9888 | truefitlab@gmail.com | truefitdentalios.com

True-Fit DENTAL IOS

Remake Scan Sheet

REMAKE CASE WORKFLOW · QUICK REFERENCE

> Every remake: notes for the lab with the specific reason and your recommended solution, plus pictures of the problem.

1 PICTURES TO TAKE FIRST

PROBLEM	PICTURE TO TAKE
Open or short margin	Crown in place
Open contact	Ideally with a measurement tool visible
Shade	Crown in the mouth, shade tab visible
Occlusion	Showing the high or light spot
Fracture	Showing the fracture

2 SCAN(S) TO TAKE

PROBLEM	SCAN(S) TO TAKE
Contacts	Adjacent teeth with the prep. **No cord needed.** Send only the scan needed.
Occlusion high or light	New **bite scan**, full and complete.
Margin or prep fit open or short margin	**Pack cord** to fully expose the margin; complete **new scan** if at all practical.*
Fracture	**Pack cord** to fully expose the margin; complete **new scan** if at all practical.*
Implant rotation, contact shift or occlusal discrepancy	New **scan body scan**: complete quadrant upper, lower and bite. Consider a new scan body.
Contour or length	Scan with the **crown in place**, showing adjacent and contralateral teeth.

*A die scan with cord packed is risky but can work if necessary.

3 SEND THE CASE WITH

- Improvements made
- Notes in the Rx
- Pictures of the problem
- The scan(s) from step 2

True-Fit Dental IOS | Trusted lab partner: Digital Dental Leaders
Eric True | **951-430-9888** | truefitlab@gmail.com | **truefitdentalios.com**

True-Fit DENTAL IOS

Remake Scan Sheet

REMAKE CASE WORKFLOW · QUICK REFERENCE

Every remake: notes for the lab with the specific reason and your recommended solution, plus pictures of the problem.

1 PICTURES TO TAKE FIRST

PROBLEM	PICTURE TO TAKE
Open or short margin	Crown in place
Open contact	Ideally with a measurement tool visible
Shade	Crown in the mouth, shade tab visible
Occlusion	Showing the high or light spot
Fracture	Showing the fracture

2 SCAN(S) TO TAKE

PROBLEM	SCAN(S) TO TAKE
Contacts	Adjacent teeth with the prep. **No cord needed.** Send only the scan needed.
Occlusion high or light	New **bite scan**, full and complete.
Margin or prep fit open or short margin	**Pack cord** to fully expose the margin; complete **new scan** if at all practical.*
Fracture	**Pack cord** to fully expose the margin; complete **new scan** if at all practical.*
Implant rotation, contact shift or occlusal discrepancy	New **scan body scan**: complete quadrant upper, lower and bite. Consider a new scan body.
Contour or length	Scan with the **crown in place**, showing adjacent and contralateral teeth.

*A die scan with cord packed is risky but can work if necessary.

3 SEND THE CASE WITH

- Improvements made
- Notes in the Rx
- Pictures of the problem
- The scan(s) from step 2

True-Fit Dental IOS | Trusted lab partner: Digital Dental Leaders
Eric True | **951-430-9888** | truefitlab@gmail.com | **truefitdentalios.com**